



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **Thermocouple boulons, acier inoxydable en fibre de verre sur cble tress - JIS Type K - Pas de filetage M12x1mm - Taille du boulon 13mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Les capteurs boulon sont généralement vissés dans la surface, ce qui augmente normalement la vitesse de réponse et la sensibilité de la mesure. Convient pour une utilisation jusqu'à +250°C.

### Specifications

## Specifications

<b>Product Code</b>	XE-0594-001
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Bolts have grounded junctions
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre stainless steel over braided cable
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	250°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Bolt thread/pitch</b>	M12 x 1mm
<b>Colour</b>	Positive: Red / Negative: White
<b>Length</b>	13mm