

**Thermocouple mis la terre fabriqué par ANSI dans un tube en acier inoxydable, câble surtress en acier inoxydable en fibre de verre - type K - 3,18 mm de diamètre - 40 mm de longueur**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

|                              |                                              |
|------------------------------|----------------------------------------------|
| <b>Product Code</b>          | XE-0634-001                                  |
| <b>General Description</b>   | 316 stainless steel tube in various sizes    |
| <b>Thermocouple Type</b>     | K                                            |
| <b>Thermocouple Junction</b> | Grounded at tip for fast response            |
| <b>Cable Length</b>          | 2M                                           |
| <b>Cable Type</b>            | Glassfibre stainless steel overbraided cable |
| <b>Cable Termination</b>     | Bare tails                                   |
| <b>Max. Temperature</b>      | 350°C                                        |
| <b>Min. Temperature</b>      | -60°C                                        |
| <b>Standards Met</b>         | ANSI MC96.1                                  |
| <b>Series</b>                | Fabricated Tube Thermocouple                 |