

**Thermocouple lames CEI, cble isol en fibre de verre avec sur-tresse en acier inoxydable - Type J avec terminaison de fiche standard - Longueur 5 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

|                              |                                                          |
|------------------------------|----------------------------------------------------------|
| <b>Product Code</b>          | XE-0716-001                                              |
| <b>General Description</b>   | Thermocouple mounted on a stainless steel shim           |
| <b>Thermocouple Type</b>     | J                                                        |
| <b>Thermocouple Junction</b> | Junction grounded at tip for fast response               |
| <b>Cable Length</b>          | 5M                                                       |
| <b>Cable Type</b>            | Glassfibre insulated lead with stainless steel overbraid |
| <b>Core / Strands</b>        | 7/0.2mm                                                  |
| <b>Termination Type</b>      | Standard Plug (with cable clamps)                        |
| <b>Max. Temperature</b>      | +350°C                                                   |
| <b>Min. Temperature</b>      | -60°C                                                    |
| <b>Dimensions</b>            | 25x13mm shim (0.2mm thick)                               |