



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

**Thermocouple de surface de tuyau ANSI,
surtresse en acier inoxydable en fibre de verre -
Type K pour sadapter un diamtre de tuyau de 30
45 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0656-001
General Description	Stainless steel pipe clamp sensor with adjustable diameter
Thermocouple Type	K
Thermocouple Junction	Grounded tip for fast response
Cable Length	2M
Cable Type	Glassfibre stainless steel braided extension cable
Core / Strands	7/0.2mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Standards Met	ANSI MC96.1
To Fit Pipe Diameters	30-45mm
Series	Pipe Surface Thermocouple