

**ANSI Type K 4.5kg Pull Power (Horseshoe)
Magnet Thermocouple, PFA Insulated Cable avec
tresse en acier inoxydable se terminant en
miniature ou standard**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	1m or 2m
Cable Type	PFA Teflon® insulated with stainless steel over-braid
Core / Strands	7/0.2mm
Termination Type	Miniature or Standard Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1