

**Sonde thermocouple usage gnral ANSI, cble isol
en fibre de verre avec sur-tresse en acier
inoxydable - Type K se terminant par une fiche
miniature - Longueur de 1 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0762-001
General Description	Rugged thermocouple with stainless steel sheath
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	1M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	ANSI MC96.1
Series	Fabricated Tube Thermocouple