

Thermocouple de buse ANSI, cble isol en fibre de verre avec sur-tresse en acier inoxydable - Type K se terminant en queues nues - Longueur 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0870-001
General Description	Thermocouple housed in a stainless steel threaded bolt
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response, M6x1mm thread
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare tails
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	ANSI MC96.1
Series	Nozzle Thermocouple