

Thermocouple rondelle ANSI, cble isol en fibre de verre avec sur-tresse en acier inoxydable se terminant par des queues nues - taille du trou de la rondelle de 10,0 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0933-001
General Description	Thermocouple housed in a washer eyelet
Thermocouple Type	K
Thermocouple Junction	Grounded tip for quick response
Cable Length	2m
Cable Type	Glassfibre insulated, stainless steel overbraided
Core / Strands	7/0.2mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Standards Met	ANSI
Washer	10.0mm clearance hole
Series	Washer Thermocouple