

Thermocouple de rondelle ANSI, cble isol en fibre de verre avec sur-tresse en acier inoxydable se terminant par une prise standard - taille du trou de la rondelle de 3,5 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0928-001
General Description	Thermocouple housed in a washer eyelet
Thermocouple Type	K
Thermocouple Junction	Grounded tip for quick response
Cable Length	2m
Cable Type	Glassfibre insulated, stainless steel overbraided
Core / Strands	7/0.2mm
Termination Type	Standard plug termination with cable clamps
Max. Temperature	+350°C
Standards Met	ANSI
Washer	3.5mm clearance hole
Series	Washer Thermocouple