

**Sonde thermocouple usage gnral CEI, cble isol
en fibre de verre avec sur-tresse en acier
inoxydable - Type J - Terminaison queue nue -
2M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-3400-001
General Description	Rugged thermocouple with stainless steel sheath
Thermocouple Type	J
Thermocouple Junction	Grounded hot junction for fast response
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare tails
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	IEC
Series	Fabricated Tube Thermocouple