




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple jonction expose IEC double torsad isol PFA avec queues nues - Type K - 1/0,315 mm - Longueur de cble de 1 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple de type K, J ou T. Rponse rapide, jonction expose soude disponible dans la spcification de cble de classe 1 conformment la norme BS EN 60584-3:2008 avec terminaison queue nue.

Ideal pour les applications de test et de developpement

Spfications

Specifications

| | |
|----------------------------|---|
| Product Code | XE-7119-001 |
| General Description | Fast response, welded exposed junction |
| Sensor Type | Exposed Junction |
| Thermocouple Type | K IEC |
| Cable Length | 1M |
| Cable Type | Teflon® PFA insulated twin-twisted lead |
| Core / Strands | 1/0.315mm |
| Cable Termination | Bare tails |
| Max. Temperature | 250°C |
| Min. Temperature | -75°C |
| Tolerance | Class 1 |
| Standards Met | IEC-584-3 |