




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple jonction expose IEC double torsad isol en PFA avec queues nues - Type J - 1/0,2 mm - Longueur de cble de 1 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple de type K, J ou T. Rponse rapide, jonction expose soude disponible dans la spcification de cble de classe 1 conformment la norme BS EN 60584-3:2008 avec terminaison queue nue.

Ideal pour les applications de test et de developpement

Spcifications

Specifications

Product Code	XE-7134-001
General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	J IEC
Cable Length	10M
Cable Type	Teflon® PFA insulated twin-twisted lead
Core / Strands	7/0.2mm
Cable Termination	Bare tails
Max. Temperature	250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3
Series	Exposed Junctions