


**Thermocouple jonction expose IEC double
torsad isol PFA avec queues nues - Type J - 1/0,5
mm - Longueur de cble de 10 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple de type K, J ou T. Rponse rapide, jonction expose soude disponible dans la spcification de cble de classe 1 conformement la norme BS EN 60584-3:2008 avec terminaison queue nue.

Idal pour les applications de test et de developpement

Spifications

Specifications

Product Code	XE-7138-001
General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	J IEC
Cable Length	10M
Cable Type	Teflon® PFA insulated twin-twisted lead
Core / Strands	1/0.5mm
Cable Termination	Bare tails
Max. Temperature	250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3