

# Thermocouple jonction expose IEC double torsion isol en PTFE avec queues nues - Type T - Longueur de cble de 3 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple de type K, J ou T. Rponse rapide, jonction expose soude disponible dans la spcification de cble de classe 1 conformment la norme BS EN 60584-3:2008 avec terminaison queue nue.

Idal pour les applications de test et de dveloppement.

Spfications

## Specifications

<b>Product Code</b>	XE-7112-001
<b>General Description</b>	Fast response, welded exposed Junction
<b>Sensor Type</b>	Exposed Junction
<b>Cable Length</b>	3M
<b>Core / Strands</b>	1/0.2mm
<b>Cable Termination</b>	Bare tails
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-75°C
<b>Tolerance</b>	Class 1
<b>Standards Met</b>	IEC-584-3