

**Sondes de temprature thermocouple visser CEI -  
Type K - 1/4" NPT - 150 mm de longueur x 3 mm  
de diamtre**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Capteurs de temperature visser de type K avec une plage de mesure de -75 +250C et un fil conducteur paire plate PFA de 1 metre.

Pour un montage direct dans des cuves, des tuyaux, des rservoirs ou des conduits pour la mesure de la temperature des liquides ou des gaz.

Applications typiques

Contrle des processus industriels

Tuyaux ou conduits

Navires et chars

### Spfications

## Specifications

<b>Product Code</b>	XE-7350-001
<b>General Description</b>	Screw in Type K thermocouple
<b>Thermocouple Type</b>	K IEC
<b>Cable Length</b>	1M
<b>Cable Type</b>	PFA
<b>Termination Type</b>	Stripped Wires
<b>Sheath Length (mm)</b>	150mm
<b>Sheath Diameter (mm)</b>	3mm
<b>Process Connection options</b>	1/4"NPT
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	IEC
<b>Number of Wires</b>	2