

Panneau aluminium pour connecteur thermocouple avec prises standard ANSI Type T - 20 voies



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Pre-assembled Standard Thermocouple Connector Panels

A pre-assembled standard horizontal panel assembly designed for use with thermocouples or leads with standard plugs. Robust and perfect for mounting in control rooms, test stands, laboratory applications, control rooms, and more. They support calibration codes including K, J, T with quick wire sockets and can accommodate 1-24 way circuits, with a maximum temperature range up to +120°C.

Supplied pre-fitted with

[Standard Fascia Sockets](#)

20 way circuit

Type T ANSI calibration

Maximum temperature 120°C

Typical Applications

OEMs

Instrumentation panel builders

Control rooms

Test Rigs and Skids

Specifications

Product Code	XE-7528-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel.
Thermocouple Type	T ANSI
Number of Ways	20 Way
Connector Size	Standard
Connector Type	Fascia Mounting Socket
Colour	Blue
Dimensions	H111mm x W221mm



Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS