



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Thermocouples tube fabriqués CEI avec joint
sertir et fil conducteur - Type K - Diamètre de la
gaine de 4,5 mm - Longueur de la gaine de 250
mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockists of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Applications typiques

Contrôle des processus industriels

Tuyaux ou conduits

Fours et chauffage

Usage général

Specifications

Specifications

Product Code	XE-7907-001
General Description	Fabricated IEC Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead
Thermocouple Type	K
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	250mm
Sheath Diameter (mm)	4.5mm
Max. Temperature	250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC
Number of Wires	2
Series	Fabricated Tube Thermocouple