



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Thermocouples tube fabriqus CEI avec joint
sertir et fil conducteur - Type J - Diamtre de la
gaine de 3 mm - Longueur de la gaine de 150 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Applications typiques
Contrôle des processus industriels
Tuyaux ou conduits
Fours et chauffage
Usage général

Specifications**Specifications**

Product Code	XE-7914-001
General Description	Fabricated IEC Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead
Thermocouple Type	J
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	150mm
Sheath Diameter (mm)	3mm
Max. Temperature	250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC
Number of Wires	2
Series	Fabricated Tube Thermocouple