

## **Thermocouples tubulaires fabriqués ANSI avec joint sertir et fil conducteur - Type K - 4,5 mm x 100 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Applications typiques**

Contrôle des processus industriels

Tuyaux ou conduits

Fours et chauffage

Usage général

**Spécifications****Specifications**

|                      |                                                                                        |
|----------------------|----------------------------------------------------------------------------------------|
| Product Code         | XE-7941-001                                                                            |
| General Description  | Fabricated Type K Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead |
| Thermocouple Type    | K                                                                                      |
| Cable Length         | 1M                                                                                     |
| Cable Type           | PFA                                                                                    |
| Termination Type     | Stripped Wires                                                                         |
| Sheath Length (mm)   | 100                                                                                    |
| Sheath Diameter (mm) | 4.5                                                                                    |
| Max. Temperature     | 250°C                                                                                  |
| Min. Temperature     | 75°C                                                                                   |
| Tolerance            | Class 1                                                                                |
| Standards Met        | ANSI MC96.1                                                                            |
| Number of Wires      | 2                                                                                      |
| Series               | Fabricated Tube Thermocouple                                                           |