



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Thermocouples tubulaires fabriqués ANSI avec joint d'étanchéité sertir et fil conducteur - Type K - 6,0 mm x 250 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockists of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Applications typiques

Contrôle des processus industriels

Tuyaux ou conduits

Fours et chauffage

Usage général

Spécifications**Specifications**

Product Code	XE-7947-001
General Description	Fabricated Type K Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead
Thermocouple Type	K
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	250
Sheath Diameter (mm)	6.0
Max. Temperature	250°C
Min. Temperature	75°C
Tolerance	Class 1
Standards Met	ANSI MC96.1
Number of Wires	2
Series	Fabricated Tube Thermocouple