

**Thermocouple IEC isol jonction fine jauge  
expose (conducteurs 0,076 mm) - Type K - 300  
mm de longueur**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple type K. Rponse rapide, jonction expose soude disponible avec queues nues. Idal pour les applications de test et dveloppement.

### Caractristiques techniques

## Specifications

<b>Product Code</b>	XE-0453-001
<b>General Description</b>	Fast response, welded tip exposed junction
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	K IEC
<b>Cable Length</b>	300mm
<b>Cable Type</b>	Teflon® PFA insulated wires
<b>Core / Strands</b>	0.003" (0.076mm) conductors
<b>Termination Type</b>	Bare Tails
<b>Max. Temperature</b>	260°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	IEC-584-3
<b>Series</b>	Exposed Junctions