




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple jonction expose isol ANSI PFA avec queues nues - Type K - Taille du brin 1/0,2 mm - Longueur 300 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0427-001
General Description	Fast response, welded exposed junction
Sensor Type	K
Thermocouple Type	Exposed Junction
Cable Length	300mm
Cable Type	Teflon® PFA insulated twin-twisted lead
Core / Strands	1/0.2mm
Termination Type	Bare Tails
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	ANSI MC96.1
Series	Exposed Junctions