



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple jonction expose IEC double torsad isol PFA avec fiche miniature - Type T - Taille du brin 7/0,2 - Longueur 300 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Jonction apparente soude réponse rapide, disponible dans la spécification de câble de classe 1 conformément à la norme BS EN 60584-3:2008 avec une terminaison de fiche miniature.

Ideal pour les applications de test et de développement.

Specifications

Specifications

Product Code	XE-0433-001
General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	T IEC
Cable Length	300mm
Cable Type	Teflon® PFA insulated twin-twisted lead
Core / Strands	7/0.2mm
Termination Type	Miniature flat pin plug colour coded in accordance IEC 584
Max. Temperature	250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3
Connector Type	Plug
Connector Size	Miniature
Series	Exposed Junctions