

**ANSI thermocouple pointe soude expose double  
torsade en fibre de verre haute temprature de  
type K avec terminaison miniature de bouchon -  
longueur 25M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractristiques techniques**
**Specifications**

<b>Product Code</b>	XE-5646-001
<b>General Description</b>	Fast response, welded exposed junction. Good temperature resistance (not suitable in fluids)
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	25M
<b>Cable Type</b>	high temperature glassfibre insulated twisted cable
<b>Core / Strands</b>	1/0.711mm conductors (21awg)
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	700°C
<b>Min. Temperature</b>	0°C
<b>Tolerance</b>	Tolerance to ANSI MC96.1 Class 2
<b>Standards Met</b>	Colour code to ANSI MC96.1
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Yellow
<b>Series</b>	Exposed Junctions