

Sonde viande sans fil RFX



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Suivez la température de la viande partir de votre appareil

Utilisez jusqu 50 sondes avec une seule passerelle

Connexion sans fil fiable

50 heures d'autonomie avec une charge de 10 minutes

Entièrement étanche selon IP69K

Certificat d'étalonnage traçable gratuit

Doit être couplé avec la
passerelle RFX *ou le*
kit de démarrage RFX

[passerelle RFX](#)
[kit de démarrage RFX](#)

Specifications

Specifications

Product Code	XE-4201-001
Model	RFX Wireless Meat Probe
Type	Thermistor, 4 Sensors
Accuracy	+/-0.5°C (-10 to 100°C)
Battery	1.5 V AAA with charging block - 175+ full charges per AAA battery. 0-80% in 10 minutes; 0-100% in 20 minutes
warranty	6 Months
Output	Wireless communication: 433Mhz radio frequency
Temperature Range	Internal: -18 to 100°C / External: -18 to 538°C
Standards Met	IP69K
Measuring range	Battery powered RFX GATEWAY: 650m (2,132 ft.), AC Powered RFX GATEWAY: 385m (1,263 ft.) (RFX MEAT to RFX GATEWAY)