



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Sondes de temperature vis ANSI - Type K - 1/4" BSPP - 200 mm de longueur x 3 mm de diamtre



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Vissez des capteurs de temprature de type K avec une plage de mesure de -75 +250C et un fil de fil paires plates PFA de 1 mtre.

Pour le montage direct dans des rcipients, tuyaux, rservoirs ou conduits pour mesurer la temprature liquide ou gazeuse.

Applications typiques

Contrle des procds industriels

Tuyaux ou conduits

Navires et Rservoirs

Caractristiques techniques

Specifications

Product Code	XE-7367-001
General Description	Screw in Type K thermocouple
Thermocouple Type	K ANSI
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	200mm
Sheath Diameter (mm)	3mm
Process Connection options	1/4"BSPP
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	ANSI
Number of Wires	2