

**Sondes de temprature thermocouple visses ANSI  
- Type K - 1/2" BSPP - 200 mm de longueur x 6 mm  
de diamtre**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Vissez des capteurs de temprature de type K avec une plage de mesure de -75 +250C et un fil de fil paires plates PFA de 1 mtre.

Pour le montage direct dans des rcipients, tuyaux, rservoirs ou conduits pour mesurer la temprature liquide ou gazeuse.

Applications typiques

Contrle des procds industriels

Tuyaux ou conduits

Navires et Rservoirs

### Caractristiques techniques

## Specifications

<b>Product Code</b>	XE-7379-001
<b>General Description</b>	Screw in Type K thermocouple
<b>Thermocouple Type</b>	K ANSI
<b>Cable Length</b>	1M
<b>Cable Type</b>	PFA
<b>Termination Type</b>	Stripped Wires
<b>Sheath Length (mm)</b>	200mm
<b>Sheath Diameter (mm)</b>	6mm
<b>Process Connection options</b>	1/2"BSPP
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	ANSI
<b>Number of Wires</b>	2