


# Connecteur miniature rectangulaire pour thermocouple douille de fascia AM-T-FF Type T ANSI



**Miniature Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Connecteur miniature rectangulaire pour thermocouple Douille fascia AM-T-FF Type T ANSI**

Clip en nylon pour une fixation rapide et simple sur le panneau. Toutes les prises de montage de la faade sont polarises au moyen de broches de taille ingale

**Caractristiques techniques****Specifications**

<b>Product Code</b>	XE-1227-001
<b>General Description</b>	All fascia mounting sockets are polarised by means of unequally sized pins
<b>Thermocouple Type</b>	T ANSI
<b>Connector Type</b>	Fascia mounting socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Blue
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS