




Cble / fil thermocouple paires plates isolés en PVC Type T (IEC) - 13/0,2 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractristiques techniques
Specifications

General Description	Good general purpose insulation for 'light' environments. Waterproof and very flexible.
Cable Length	10, 25, 50, 100 & 200m
Cable Type	Flat pair construction
Core / Strands	13/0.2mm
Conductor Cross Section 13/0.2 (mm²)	0.408
Outer Colour	Brown
+ Colour	Brwon
- Colour	White
Drain Wire	No
Insulation Material	PVC Insulated
Loop Resistance 13/0.2 (Ohms per combined metre)	1.2
Max. Temperature	105°C
Min. Temperature	-10°C
+ Conductor	Copper
- Conductor	Constantan
Screen	No
Standards Met	Colour code (cores & jacket) to IEC-584-3



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termination Type	T
Tolerance	Class 2