


Prise de connecteur thermocouple pour montage de PCB AM-T-PCB Type T ANSI - broches plaques tain



Miniature PCB Socket
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Prise de connecteur thermocouple pour montage de PCB AM-T-PCB-TIN Type T ANSI

Spécialement conçus pour les circuits imprimés, broches en tain massif pour le montage direct sur carte et la fourniture d'un capteur CJC

Caractéristiques techniques

Specifications

Product Code	XE-1230-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	T ANSI
Connector Size	Miniature
Connector Type	PCB socket
Colour	Blue
Max. Temperature	220°C
Pins	Tin Plated
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS