



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple ISOL Minral JIS avec bouchon standard - Type K - 1,5 mm x 500 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Vous avez besoin d'aligner votre capteur de température ?

Traable selon les normes d'alignement ISO17025 (UKAS), Labfacility propose le choix entre un alignement traable en 3 ou 5 points via notre installation d'alignement interne.

Veuillez noter que si vous ajoutez une option de calibration votre commande, vous serez contacté pour les informations de livraison.

Caractéristiques techniques**Specifications**

Product Code	XE-2916-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction insulated from sheath (offers protection against spurious electrical signals)
Probe Diameter	1.5mm
Probe Length	500mm
Sheath Material	310 Stainless Steel
Termination Type	Standard round pin plug, colour coded in accordance JIS
Max. Temperature	1100°C
Min. Temperature	-40°C
Standards Met	JIS
Connector Type	Plug



Connector Size	Standard
Colour	Blue
Series	Mineral Insulated