

Thermocouple minral isol JIS avec bouchon standard - Type K - 3,0 mm x 250 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Vous avez besoin d'alonner votre capteur de température ?

Traable selon les normes d'atlonnage ISO17025 (UKAS), Labfacility propose le choix entre un talonage traable en 3 ou 5 points via notre installation d'atlonnage interne.

Veillez noter que si vous ajoutez une option de calibration votre commande, vous serez contact pour les informations de livraison.

Caractéristiques techniques**Specifications**

Product Code	XE-2914-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction insulated from sheath (offers protection against spurious electrical signals)
Probe Diameter	3.0mm
Probe Length	250mm
Sheath Material	310 Stainless Steel
Termination Type	Standard round pin plug, colour coded in accordance JIS
Max. Temperature	1100°C
Min. Temperature	-40°C
Standards Met	JIS
Connector Type	Plug



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Connector Size	Standard
Colour	Blue
Series	Mineral Insulated