



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple ANSI pointe expose double extrmit soude torsade haute temprature ANSI avec terminaison miniature de bouchon - longueur de 1 M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractristiques techniques
Specifications

Product Code	XE-5647-001
General Description	Fast response, welded exposed junction. Good temperature resistance (not suitable in fluids).
Sensor Type	Exposed Junction
Thermocouple Type	K
Cable Length	1M
Cable Type	high temperature glassfibre insulated twisted cable
Core / Strands	1/0.711mm conductors (21awg)
Cable Termination	Miniature Plug
Max. Temperature	700°C
Min. Temperature	0°C
Tolerance	Tolerance to ANSI MC96.1 Class 2
Standards Met	Colour code to ANSI MC96.1
Connector Type	Plug
Connector Size	Miniature
Colour	Yellow
Series	Exposed Junctions