



Prise de connecteur thermocouple de montage de PCB AM-E-PCB Type E ANSI



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Prise de connecteur thermocouple de montage de PCB AM-E-PCB Type E ANSI

Spécialement conçus pour les circuits imprimés, broches pleines pour le montage direct sur carte et provision pour capteur CJC

Caractéristiques techniques**Specifications**

Product Code	XE-9039-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	E ANSI
Connector Type	PCB socket
Connector Size	Miniature
Colour	Violet
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS