


Prise standard pour connecteur de thermocouple rectangulaire fascia AS-R/S-FF Type R/S ANSI



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur de thermocouple rectangulaire standard Prise fascia AS-R/S-FF Type R/S ANSI

Clip en nylon pour une fixation rapide et simple sur le panneau. Toutes les prises de montage de la faade sont polarises au moyen de broches de taille ingale

Caractristiques techniques**Specifications**

Product Code	XE-1617-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	R/S ANSI
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Green
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS