


Connecteur miniature de thermocouple Douille en ligne JM-E-FS Type E JIS



Miniature in-line Socket
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur miniature de thermocouple Douille en ligne JM-E-FS Type E JIS

Ideal pour un usage gnral, tous les contacts sont polariss pour garantir une bonne connexion. Ces connecteurs acceptent des cbles thermocouple jusqu 4 mm de diamtre.

Caractristiques techniques**Specifications**

Product Code	XE-9232-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	E JIS
Connector Type	Socket
Connector Size	Miniature
Colour	Violet
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS