


## Fichire de connecteur thermocouple fil rapide miniature JM-T-MQ Type T JIS



**Miniature in-line Plug - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Connecteur miniature pour thermocouple, fiche rapide JM-T-MQ Type T JIS**

Connecteur miniature avec broches de prise plates. Connexion facile par injection pour une limitation rapide. Il suffit de pousser le fil et de serrer la vis.

**Caractéristiques techniques****Specifications**

<b>Product Code</b>	XE-9262-001
<b>General Description</b>	Easy jab-in connection for rapid termination.
<b>Thermocouple Type</b>	T JIS
<b>Connector Type</b>	Quick Wire Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Brown
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Flat Pins
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS