


Prise de connecteur de thermocouple fil rapide miniature JM-T-FQ Type T JIS



Miniature in-line Socket - Quick Wire
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur fil rapide miniature pour douille thermocouple JM-T-FQ Type T JIS

Connexion facile piquer, il suffit de foncer le fil et de serrer la vis. Les contacts sont polarisés pour éviter une terminaison incorrecte.

Caractéristiques techniques**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-9263-001 |
| General Description | Easy jab-in connection, just push in wire & tighten screw. |
| Thermocouple Type | T JIS |
| Connector Type | Quick Wire Socket |
| Connector Size | Miniature |
| Colour | Brown |
| Max. Temperature | 220°C |
| +Positive Contact | Copper |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |