


Prise miniature pour connecteur thermocouple fil rapide JM-R/S-FQ Type R/S JIS



Miniature in-line Socket - Quick Wire
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur fil rapide miniature Douille thermocouple JM-R/S-FQ Type R/S JIS

Connexion facile piquer, il suffit d'enfoncer le fil et de serrer la vis. Les contacts sont polarisés pour éviter une terminaison incorrecte.

Caractéristiques techniques**Specifications**

Product Code	XE-9203-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	R/S JIS
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	Black
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS