


Connecteur de thermocouple douille de fascia miniature rectangulaire JM-E-FF Type E JIS



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur miniature rectangulaire pour thermocouple Douille de fascia JM-E-FF Type E JIS

Clip en nylon pour une fixation rapide et simple sur le panneau. Toutes les prises de montage de la faade sont polarises au moyen de broches de taille ingale

Caractristiques techniques**Specifications**

Product Code	XE-9237-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	E JIS
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	Violet
Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS