




Prise de connecteur thermocouple monte sur PCB JM-B-PCB Type B JIS



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Prise de connecteur thermocouple monte sur PCB JM-B-PCB Type B JIS

Spécialement conçus pour les circuits imprimés, broches pleines pour le montage direct sur carte et provision pour capteur CJC

Caractéristiques techniques**Specifications**

Product Code	XE-9220-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	B JIS
Connector Type	PCB Socket
Connector Size	Miniature
Colour	Grey
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS