

Standard Thermocouple Duplex Thermocouple Connecteur Plug JS-K-MD Type K JIS



Standard Duplex in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard Thermocouple Duplex Thermocouple Connecteur Plug JS-K-MD Type K JIS

Corps en polyester remplis de verre accidenté. Bouchons équipés d'épingles thermocouples solides. Les contacts sont polarisés pour éviter une connexion incorrecte

Spécifications**Specifications**

| | |
|----------------------------|---|
| Product Code | XE-9286-001 |
| General Description | Rugged glass filled polyester bodies, suitable for industrial applications. |
| Thermocouple Type | K JIS |
| Connector Type | Duplex Plug |
| Connector Size | Standard |
| Colour | Blue |
| Dimensions | See Schematic Drawing (Product Literature Tab) |
| Max. Temperature | 220°C |
| Pins | Solid Round Pins |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Nickel Alloy |
| Standards Met | BS EN 50212:1997. RoHS |